

HFS MCRIF32 System Requirements

The minimum system requirements for the HFS 32-bit software

Processor (CPU)

1.5 GHz or higher – x86 or x64 compatible processor.

Memory (RAM)

1.5 gigabytes (GB) or more for 32-bit OS (Recommend 3GB or more).

4 gigabytes (GB) or more for 64-bit OS (Recommend 8GB or more).

Hard Drive

1.5 gigabytes (GB) of free space.

Writable Media

Floppy, USB drives and/or CD/DVD writable.

Interface

Keyboard and mouse or some other compatible pointing device.

Display

1024x768 or higher resolution monitor.

Hard Drive

1.5 gigabytes (GB) of free space.

Operating System

Microsoft Windows 7/8/8.1/10, or later operating systems.

Installation and Registry

MCRIF32 installation requires an account with administrative rights. During the installation of the software, the following registry keys are modified:

HKEY_LOCAL_MACHINE\SOFTWARE\Classes\CLSID

HKEY_LOCAL_MACHINE\SOFTWARE\Classes\CLSID

HKEY_LOCAL_MACHINE\SOFTWARE\Classes\Installer\Features
HKEY_LOCAL_MACHINE\SOFTWARE\Classes\Installer\Products
HKEY_LOCAL_MACHINE\SOFTWARE\Classes\Installer\UpgradeCodes
HKEY_LOCAL_MACHINE\SOFTWARE\Classes\Installer\Win32Assemblies
HKEY_LOCAL_MACHINE\SOFTWARE\Classes\Interface
HKEY_LOCAL_MACHINE\SYSTEM\ControlSet001\Enum\Root\LEGACY_MSISERVER\0000\Control
HKEY_LOCAL_MACHINE\SYSTEM\CurrentControlSet\Enum\Root\LEGACY_MSISERVER\0000\Control
User's settings are saved in HKEY_CURRENT_USER\Software\Health Financial Systems key

Others

Internet Explorer 7.0 or later, 32 bit browser only. Internet functionality requires Internet access (fees may apply).

Microsoft XML Core Services (MSXML) version 3.0 and 6.0.

Adobe Reader version 8 or higher.

Files and Folders Permissions

Software Installation Folder - The software application creates temporary files and folders in the application installation folder, therefore the User Account must have permission to read, write, create, delete and modify files and subfolders in the folder where the HFS Cost Report software application is installed. These permissions are equivalent to the "Full Control" permissions setting.

Data Files Folder - The software application creates files for the purpose of storing cost report data, therefore the User Account must have permission to read, write, create, delete and modify files in a folder you designate as the location of the user data files. The data files folder may be the same of different than the software application installation folder. These permissions are equivalent to the "Full Control" permissions setting.

Internet Connection

ECR Export and Check for Updates require Internet connection (fees may apply).

Additional

Actual requirements and product functionality may vary based on your system configuration and operating system.